



## GENOTYPING – Bridgewater Laboratory Protocol

### Preparations:

- Retrieve the samples from the CAF freezer. Put the samples in a -20°C freezer until ready for genotyping.
- Make sure all reagents and materials are present.
  - ddH<sub>2</sub>O
  - 5X Lysis Buffer A (-20°C – box “Genotyping”)
  - 10X Lysis Buffer B (-20°C – box “Genotyping”)
  - 2X HS-Red TAQ Mix (-20°C – box “Genotyping”)
  - Primers (-20°C – box “Genotyping”)
  - RedSafe (+4°C box “Genotyping”)
  - 1 kb and/or 100 bp DNA ladder(s) (+4°C box “Genotyping”)
  - Agarose (chemicals shelf)
  - Small Erlenmeyer flask (bench where electrophoresis is performed)
  - 100 ml cylinder (bench where electrophoresis is performed)
  - Pipettes and pipette tips
- Write down the specific protocol and make the calculations beforehand.
  - Which samples (the written numbers as indicated by the CAF; write down how many samples there are in total and their numbers) need to be processed.
  - Which primer to use for each sample.
  - Count the number of samples for each primer and the total number of samples and plan how to lay out the samples on the agarose gel – the gel has 20 columns and 2 rows.
  - Keep in mind to include at least one column for a ladder for each side of the gel, one column for the positive control, and one column for the negative control for each primer being genotyped for.
- **Make a schedule about when to perform each step.**
  - a) Lysis of samples: ~1h
  - b) Set up PCR: ~30min.
  - c) PCR will run for ~ 2hrs
  - d) Make the agarose gel: ~40min
  - e) Run gel electrophoresis + take picture of the gel + upload results: ~1h
    - DAY 1: a) + b) = total of ~1.5hrs → then run c) = total of ~2hrs OR overnight (proceed to *DAY 2* if protocol is carried on to the next day)
    - DAY 2: d) + e) = total of ~2hrs



### Step 1 – Lysis

- 1) Get samples out from -20°C freezer and line up on a 1.7 ml tube rack in numerical order.
- 2) Take out 5X Buffer A (white cap) + 10X Buffer B (red cap) from -20°C freezer & thaw.
- 3) Get 50ml tube with ddH<sub>2</sub>O (from *Milli-Q*).
- 4) Add **70 µl of ddH<sub>2</sub>O** to each sample.
- 5) Add **20 µl of 5X Buffer A** to each sample.
- 6) Add **10 µl of 10X Buffer B** to each sample.

\*\*\* *TIP –instead of pipetting 70ul of ddH<sub>2</sub>O, 20ul of 5X Buffer A, and 10ul of 10X Buffer B one at a time, follow this formula:*

$$70 \text{ ul ddH}_2\text{O} \times (x) = \text{\_ul}$$

$$20 \text{ ul 5X Buffer A} \times (x) = \text{\_ul}$$

$$10 \text{ ul 10X Buffer B} \times (x) = \text{\_ul}$$

Where (x) is the total number of samples. For example, for 20 samples, (x) will equal to 20 samples.

- 7) Vortex the samples and lysis buffer contents well. Centrifuge on *MBI centrifuge*.
- 8) Incubate on the *MBI Thermo-Shaker* for 5 minutes at 95°C with 1000 speed (shaking).
  - Open the tubes' caps a few times, close the caps again when using the vortex. Vortex well. Centrifuge samples on *MBI centrifuge*.
- 9) Put back on the *MBI Thermo-Shaker* at 95°C, with 1000 speed (shaking), and incubate for another 10 minutes.
  - Open the tubes' caps a few times, close the caps again when using the vortex. Vortex well. Centrifuge samples on *MBI centrifuge*.
- 10) Put the samples back on the rack and add **900 µl of ddH<sub>2</sub>O** to each sample.
- 11) Put samples in *Centrifuge 5424* and spin for 2 minutes at 12500 rcf (or rpm).
- 12) Put samples back on the rack and store in the 4°C fridge until further use.

### Step 2 – Polymerase Chain Reaction (Set up ONE PCR machine for each primer)

- 13) Get ice.
- 14) Label PCR tubes for the samples, positive control, and negative control.
- 15) Label ONE Eppendorf tube for each primer Master Mix (i.e., for each kind of primer) and put on ice.
- 16) Take out 2X HS-Red TAQ mix and the primers that will be used from the -20°C freezer. Thaw (on ice).
- 17) Make the primer master mix for each primer (on ice).



- **5 µl 2X HS-Red TAQ mix /sample** → Multiply by number of samples + 1 positive control + 1 negative control
- **3.6 µl ddH<sub>2</sub>O / sample** → Multiply by number of samples + 1 positive control + 1 negative control
- **0.4 µl Primer Mix / sample** → Multiply by number of samples + 1 positive control + 1 negative control

18) Set the labeled PCR tubes onto purple *isofreeze* box (in the bottom shelf of -20°C freezer).

19) Add **9 µl of primer master mix** to each PCR tube.

20) Add **1 µl of each sample** to each PCR tube (get a POSITIVE control from +4°C box “Genotyping”)

21) Close the lid of the PCR tubes tightly.

22) Vortex the PCR tubes.

23) Spin the PCR tubes down in the *Hoefer mini-centrifuge* to get rid of most of the bubbles and bring the contents back to the bottom of the tube.

24) Put PCR tubes into the *MBI PCR machine(s)* and start the program(s).

LEFT <i>MBI PCR MACHINE</i> (#1)	RIGHT <i>MBI PCR MACHINE</i> (#2)
<ul style="list-style-type: none"> <li>• <b>Cre</b> → ‘login’ under “Antje” → ‘enter’ → ‘program’ → “cre genotyping” → ‘start’</li> <li>• <b>β-gal (Shrm3)</b> → ‘login’ under “Antje” → ‘enter’ → ‘program’ → “shrm3/bgal” → ‘start’</li> <li>• <b>Shroom3-LoxP</b> → ‘login’ under “Kristina” → ‘enter’ → ‘program’ → “shroom3-loxp” → ‘start’</li> <li>• <b>Jarid F/R</b> → ‘login’ under “Kristina” → ‘enter’ → ‘program’ → “jarid sexing” → ‘start’</li> <li>• <b>Six2Cre-Insertion</b> → ‘login’ under “Kristina” → ‘enter’ → ‘program’ → “six2cre-insert” → ‘start’</li> </ul>	<ul style="list-style-type: none"> <li>• <b>LoxP</b> → ‘login’ under “tony” → ‘enter’ → ‘program’ → “loxp-new” → ‘start’</li> <li>• <b>Cre</b> → ‘login’ under “Felix” → ‘enter’ → ‘program’ → “genotyping” → ‘start’</li> <li>• <b>Shroom3-LoxP</b> → ‘login’ under “Kristina” → ‘enter’ → ‘program’ → “shroom3-loxp” → ‘start’</li> <li>• <b>Jarid F/R</b> → ‘login’ under “Kristina” → ‘enter’ → ‘program’ → “jarid sexing” → ‘start’</li> <li>• <b>Six2Cre-Insertion</b> → ‘login’ under “Kristina” → ‘enter’ → ‘program’ → “six2cre-insert” → ‘start’</li> </ul>

25) Wait ~2hrs OR let the PCR run overnight and continue the protocol the next day, i.e., Day 2.



### Step 3 – Gel Electrophoresis

- 26) Measure **100ml of 1X TAE buffer** and add to Erlenmeyer flask.
- 27) Weigh out **1g of agarose** and add to Erlenmeyer flask. This 1% agarose applies to *all* primers being used in the laboratory *EXCEPT* for the Jarid F/R primer.
  - IF genotyping for sex (using the Jarid F/R primers), weigh out **2g of agarose** (agarose is on the chemicals shelf above the scale) and add to Erlenmeyer flask. This 2% agarose applies to **ONLY** the Jarid F/R primer.
- 28) Microwave for 30 seconds → stop and swirl → microwave for 30 seconds → stop and swirl → microwave for 15 to 20 seconds (stop when it starts to boil) – CAUTION: Erlenmeyer flask will be HOT.
  - IF genotyping for sex (using the Jarid F/R primers), add an extra 15 to 20 seconds (stop when it starts to boil). This **ONLY** applies to the 2% agarose gel.
- 29) Add **5 µl of RedSafe** to Erlenmeyer flask and swirl until the agarose gel is homogenous.
- 30) Put the gel cast into the apparatus and make sure it is tight.
- 31) Pour in the 1g agarose and 100ml TAE buffer solution into the gel cast.
- 32) Put in the gel combs.
  - Put the thicker side of the 20-well combs down.
- 33) Wait ~ 30 minutes for the gel to solidify
- 34) Take out the combs carefully (upwards motion).
- 35) Transfer the gel cast over to the *OWL* gel chamber. Check that the TAE buffer volume is up to the ‘fill line’.
- 36) Take samples out from the *MBI PCR machine(s)*.
  - Press ‘STOP’ make sure the program has stopped.
  - Select ‘EXIT’ or ‘LOG OUT’
  - Turn the *MBI PCR machine* OFF.
- 37) Load **10 µl of 1 kb ladder** in well #1 on each side of the gel.
- 38) Load **10 µl of each sample** from each PCR tube into corresponding well.
  - When loading, be careful **NOT** to damage or poke the gel. Release the contents from the pipette tip *just above* the well.
- 39) Close the *OWL* gel chamber lid and turn on *BioRad*. Adjust the settings to ‘Run’ at 135V for 45 to 60 minutes (depending on the primers being used).
  - IF genotyping for sex (using the Jarid F/R primers), set the time to at least 60 minutes.
  - IF genotyping other primers (**NOT** the Jarid F/R primers), set the time to 50 minutes.



40) Carefully take out the gel cast with the agarose gel from the *OWL* gel chamber.

41) Take a picture of the gel.

- Put the agarose gel inside the *SYNGENE* machine (NOT including the gel cast).
- Click on the ‘**GeneSys**’ application in the computer.
- Click on ‘**Ethidium Bromide**’.
  - Only click *after* the interface on the left side of the screen (where it says gels, blots, etc.) goes from black/white to coloured.
- Modify zoom/brightness of the gel as needed.
- ‘Invert’ the picture to make the samples black and the background white/light gray).
- Save the image and print a copy → click the printer icon in the page that pops out.
- Click “X” to close program
- Take out the agarose gel and gently clean the *SYNGENE* machine with ddH<sub>2</sub>O and paper towels.

42) Clean up.

- Throw the agarose gel into biohazardous garbage (yellow bags).
- Close the *OWL* gel chamber and tightly cover with plastic wrap.
- Make sure all the reagents and samples are put back into their respective places, i.e., the fridge or the freezer.
- Wash any used materials (scoops, glassware, graduated cylinders, etc.).

43) Enter the genotyping results in OneDrive for the CAF staff.

- For samples that used Cre and Six2Cre-Insertion primers, upload results as “**Cre Pos**” or “**Cre Neg**”
- For samples that used β-gal (Shroom3), upload results as “**WT**” or “**Het**”
- For samples that used Shroom3-LoxP, upload results as “**LoxP**”
- For samples that used Jarid F/R, upload results as “**MALE**” or “**FEMALE**”

44) Send an e-mail to the CAF technicians to inform them about the uploaded genotyping results on OneDrive.

Attach the link to the OneDrive Sheet in the e-mail.

- Include the actual results in the text of the email as well, as this is some of the staff’s preference, ex., **261 → CRE POS**